






























7.03_Bottlelist water

Article no. name	material cap colour	volum e	GHS-label	containing preservation filling instruction	No. Bottles, Analyseparameter
 A004 Neutral	PE, white red	500 ml			1 I ⁻ , Br ⁻ , BrO ₃ ⁻ , Cl ⁻ , Cr(VI), Extinction/color, F ⁻ , NO ₂ ⁻ , NO ₃ ⁻ , o-PO ₄ ³⁻ , Silicate, SO ₄ ²⁻ , TOC/DOC, turbidity, KMnO ₄ 1 PFC, Glyphosate/AMPA, EDTA/NTA, Dry residue 1 Glowbacteria 1 BOD 2 solid material, settable solids 1 daphnia test 1 fishegg test 1 green alga test 2 alpha-activity, total 2 beta-activity, total 4 waterlens test 2 Tritium for age determination 1 Tritium for drinking water ordinance
 A703 Neutral	PE, white red	1 L			1 filterable substances 2 settleable solids
 A400 Organics	glass, green black	500 ml		fill only 90%	1 EOX, odor, anionic, cationic and non ionic detergents 1 PAH 1 PCB 2 TPH (GC) 1 Epichlorhydrin 3 Phthalates 3 Pesticides / plant protection agents 1 Pharmaceuticals / Sweeteners 2 Organotin compounds
 A002 Mibio	PE, white blue	250 ml		Thiosulphate don't rinse, fill completely – don't overrun	1 E.coli, coliforms, bacterial count 20°C and 36°C + x (x = Enterokoccus or Cl.perfringens or Ps.aeruginosa) 1 For each additional parameter
 A207 Mibio (pools)	PE, white blue	250 ml		Thiosulphate don't rinse, fill completely – don't overrun	1 Scooped sample see A002
 A001 Legionella	PE, white blue	125 ml		Thiosulphate don't rinse, fill completely – don't overrun	1 Legionella
 A203 CO₂	PET, clear white	250 ml		blasenfrei bis zum Rand	1 Bc 1 Ac., HCO ₃ ⁻ , CO ₂ , CO ₃ ²⁻ , pH, conductivity 1 EDTA/NTA
 A200 AOX	glass, brown green	250 ml		HNO ₃ don't rinse	1 AOX 4 Öl and Fat (IR)

For further information, please contact your responsible CRM in the laboratory. Preservation according to EN ISO 5667-3.



7.03_Bottlelist water

Article no. name	material cap colour	volum e	GHS-label	containing preservation filling instruction	No. Bottles, Analyseparameter
 A208 COD, N, P	PE, white white	250 ml		H ₂ SO ₄ don't rinse	1 COD, TNb, TKN, P, NH ₄ ,
 A102 Metals	PE, white black	100 ml		HNO ₃ don't rinse	1 B, metals, silicates (Attention! ground water: filtered filling!)
 A103 VOC	glass, clear silver	20 ml		H ₂ SO ₄ do not rinse out, fill bottles to the brim	2 VOC fully filled (BTX/LHKW for seepage, ground, surface and waste water) 1 VOC half filled (BTX/LHKW for seepage, ground, surface and waste water) 2 Methane, ethane, ethylene
 A101 VOC (pools)	glass, light brown silver	20 ml		Thiosulphate don't rinse, fill all to the brim, for bathing water add 3 drops of KHSO ₄ solution	2 VOC fully filled (BTX/LHKW for seepage, ground, surface and waste water)
 A201 Rn	glass, brown black	250 ml		bubblefree to the brim	1 Radon 1 Acrylamid 1 Glycols
 A104 Chlorate, Chlorite	PP, clear white	30 ml		NaOH don't rinse	1 Chlorite, Chlorate
 A211 Sulphides	PE, white green	250 ml	 	ZnAc don't rinse adjust to pH 8,5- 9,0 with NaOH	1 Sulphides
 A105 Sulphites	PE, white white	60 ml	 	EDTA don't rinse	1 Sulphites
 A106 Phenolindex	glass, brown white	60 ml		H ₂ SO ₄ don't rinse	1 Phenolindex
 A401 Phenoles	glass, green white	500 ml		H ₂ SO ₄ don't rinse	1 Phenoles

7.03_Bottlelist water

Article no. name	material cap colour	volum e	GHS-label	containing preservation filling instruction	No. Bottles, Analyseparameter
------------------	---------------------	---------	-----------	---	-------------------------------

For further information, please contact your responsible CRM in the laboratory. Preservation according to EN ISO 5667-3.

	A107 Fe II, Mn II Hg	glass, clear black	100 ml		HCl don't rinse bubblefree to the brim	1 Fe (II) 1 Mn (II) 1 Hg (Attention! ground water: filtered filling!)
	A108 Cyanide	PE, white green	60 ml		NaOH don't rinse	1 Cyanide
	A009 Chlorophyll	PE, black black	500 ml			2 Chlorophyll
	A109 Oxygen	PET, clear white	125 ml		bubblefree to the brim, add solution 1+2	1 Oxygen
	A900 O ₂ fixing solutions	glass, brown white-red	100 ml		MnCl ₂ (Solution 1) alkaline Iodide- Azide-Solution (Solution 2)	Solutions for O ₂ -fixation in A109 1 Solution 1 1 Solution 2
	A204 CO ₂ marble test	PET, clear black	250 ml		Marble don't rinse	1 CO ₂ (chalk dissolves)
	A704 Metalle UBA	PE, white, black	1 L		HNO ₃ nicht ausspülen	1 Metals by UBA
	A800 Radio	PE, white white	5000 ml			1 Radioactivity

For further information, please contact your responsible CRM in the laboratory. Preservation according to EN ISO 5667-3.

7.03_Bottlelist water