



























7.03_Bottlelist water and drinking water

Article no. name	material cap colour	volume	GHS-label	containing preservation filling instruction	No. Bottles - parameter
 A004 Neutral	PE, white red	500 ml			1 I, Br, BrO ₃ ⁻ , Cl ⁻ , Cr(VI), Extinction/color, F ⁻ , NO ₂ ⁻ , NO ₃ ⁻ , o-PO ₄ ³⁻ , SO ₄ ²⁻ , TOC/DOC, turbidity, KMnO ₄ 1 PFC, Glyphosate/AMPA, Dry residue 1 Glowbacteria 2 BOD 1 daphnia test 1 fishegg test 1 green algea test 2 alpha-activity, total 2 beta-activity, total 4 waterlens test 2 Tritium for age determination 1 Tritium for drinking water ordinance
 A703 Neutral	PE, white red	1 L			1 solid material 2 settable solids
 A400 Organics	glass, green black	500 ml		fill only 90%	1 EOX, odor 1 anionic, cationic and non ionic detergents 1 PAH, 1 PCB, 1 TPH (GC) 3 Phthalates 3 Pesticides / plant protection agents 1 Drugs / Sweeteners 2 Organotin
 A002 Mibio	PE, white red	250 ml		Thiosulphate don't rinse, fill completely – don't overrun	1 E.coli, coliforms, bacterial count 20°C and 36°C + x (x = Enterokoccus or Cl.perfringens or Ps.aeruginosa) 1 For each additional parameter
 A207 Mibio (pools)	PE, white blue	250 ml		Thiosulphate don't rinse, fill completely – don't overrun	Scooped sample see A002
 A001 Legionella	PE, white blue	125 ml		Thiosulphate don't rinse, fill completely – don't overrun	1 Legionella
 A203 CO ₂	PET, clear white	250 ml		bubblefree to the brim	1 Bc, Ac, HCO ₃ ⁻ , CO ₂ , CO ₃ ²⁻ , pH, conductivity, EDTA/NTA
 A200 AOX	glass, brown green	250 ml		HNO ₃ don't rinse	1 AOX 4 Oil and Fat (IR)

For further information please ask your responsible CRM. All preservations according to EN ISO 5667-3 and additional local criteria

Erstellt	Tanja Plaumann	24.10.2019
Geprüft	Sibylle Steinbeiß	24.10.2019
Freigegeben	Tanja Plaumann	24.10.2019











7.03_Bottlelist water and drinking water

	A208 COD, N, P	PE, white white	250 ml		H ₂ SO ₄ don't rinse	1	COD, TNb, TKN, P, NH ₄ ,
	A102 Metals	PE, white black	100 ml		HNO ₃ don't rinse	1	B, metals, silicate
	A103 VOC	glass, clear silver	20 ml		H ₂ SO ₄ don't rinse Fill all to the brim	2 2	VOC Methane, Ethane Ethene
	A101 VOC (pools)	glass, light brown silver	20 ml		Thiosulphate don't rinse, fill all to the brim, for bathing water add 3 drops of KHSO ₄ solution	2	VOC, (BTX/LHKW for drinking and bathing water), THM
	A201 Rn	glass, brown black	250 ml		bubblefree to the brim	1 1 1	Radon Epichlorohydrin Glykole
	A104 Chlorate, Chlorite	PP, clear white	30 ml		NaOH don't rinse	1	Chlorite, Chlorate
	A211 Sulfide	PE, white green	250 ml		ZnAc+NaOH don't rinse	1	Sulfide
	A105 Sulfite	PE, white white	60 ml		EDTA don't rinse	1	Sulfit
	A106 Phenolindex	glass, brown white	60 ml		H ₂ SO ₄ don't rinse	1	Phenolindex
	A401 Phenoles	glass, green white	500 ml		H ₂ SO ₄ don't rinse	1	Phenoles
	A107 Fe II, Mn II	glass, clear black	100 ml		HCl don't rinse bubblefree to the brim	1 1 1	Fe (II), Mn (II) Hg (groundwater : filter sample)

For further information please ask your responsible CRM. All preservations according to EN ISO 5667-3 and additional local criteria

Erstellt	Tanja Plaumann	24.10.2019
Geprüft	Sibylle Steinbeiß	24.10.2019
Freigegeben	Tanja Plaumann	24.10.2019

7.03_Bottlelist water and drinking water

	A108 Cyanide	PE, white green	60 ml		NaOH don't rinse	1	Cyanide
	A009 Chlorophyll	PE, black black	500 ml			2	Chlorophyll
	A109 Oxygen	PET, clear white	125 ml		bubblefree to the brim, add solution 1+2	1	Oxygen
	A900 O ₂ fixing solutions	Glas, brown white- red	100 ml		MnCl ₂ (Solution 1) alkaline Iodide-Azide- Solution (Solution 2)	1 1	Solutions for O ₂ -fixation in A109 Solution 1 Solution 2
	A204 CO ₂ marble test	PET, clear black	250 ml		Marble don't rinse	1	CO ₂
	A704 Metalle UBA	PE, white blue	1 L		HNO ₃ don't rinse	1	Metals by UBA
	A800 Radio	PE, white white	5000 ml			1	Radioactivity

For further information please ask your responsible CRM. All preservations according to EN ISO 5667-3 and additional local criteria

Erstellt	Tanja Plaumann	24.10.2019
Geprüft	Sibylle Steinbeiß	24.10.2019
Freigegeben	Tanja Plaumann	24.10.2019