

material order form water, DIN EN ISO/IEC 17025:2018, chapter 7.04

Responsible for the MF: ALBR CRM 60 – Head of department

Document user: extern (agrolab.com), ALBR CRM 10 / 60 (with subunits)

Related QM-document: MA-00667-DE, MF-02494-DE, MF-02493-DE, MF-02495-DE

Please send this form by fax or email to: **+49 8765 93996-89** or versand.bruckberg@agrolab.de

Please fill out one order form each for different sets.

Customer name:

Project:

Customer ID.:

Serviceteam AGROLAB (if known):

Contact person:

Street

ZIP City

REQUIERED FIELDS *

*Date of order

* Phone number for quiries

* Required until

Number of sets:	Sets per box:
------------------------	----------------------

Number of bottles per set	Article no	Name	Accessories (please specify quantity per set)
	A004	Neutral	
	A700	Organics 1l (TPH, PAH, PCB, pesticides)	for A120/107 Disposable filter for A120/107 yellow
	A205	Organics (alcohols, glycols)	for A060/A115 Disposable filter for A060/A115 white
	A400	Organics	for A060/A115 Disposable syringes 20ml
	A401	Organics H ₂ SO ₄	for A060/A115 Disposable syringes 60ml
	A103	VOC	A997/ A for A703 with addition B for A700 box(es) with cooling elements
	A112	VOC-P&T	A997/ B
	A113	VOC 100ml	A998 Box(es) with cooling elements
	A102	Metals	UPS Return Label
	A107	Hg	
	A119	Metals, HCl (Fe II, Mn II, Sn)	
	A120	Metals, filtered	
	A201	Rn	
	A203	CO ₂	
	A204	CO ₂ marble test	
	A115	Neutral, filtered	
	A703	Neutral 1l	
	A060	P,N, DOC filtered	
	A106	Phenolindex	
	A200	AOX	
	A208	COD, N,P	
	A104	Chlorate, Chlorite	
	A114	Cyanide	
	A210	Cyanide (waste water)	
	A211	Sulfide	
	A105	Sulfite	
	A002	Mibio for Eching	
	A009	Chlorophyll	
	A109	Oxygen	
	A 900	Fix 1+ Fix 2 for O ₂ (in add. to A109)	

Notes:

In the AGROLAB **bottle set list** you can find all the information regarding our bottles and the associated parameters. Take a look at our **handbook** for the "correct use of the AGROLAB bottle set and transport system". Here, you can find information on the correct use of our transport system for safe and undamaged transport of your samples to our laboratory. You can find all documents in the download area on our [website](#).
 For the delivery of samples to AGROLAB, please only use the bottles that we specify for this purpose, since we can only ensure the quality of sampling, delivery and subsequent analysis in this way. **Please do not cover the barcode already existing on the bottles.** The bottles remain the property of AGROLAB.

To be completed by the laboratory:

Shipping by:

Remarks customer

Version created	27.06.2023	Anja Grubert
Reviewed	27.06.2023	Mandy Erdmann-Schiessling
Approved	07.07.2023	Katja Opitz